*Speed* Biogram™

English	.....	1
Français	.....	7
Español	.....	13
Português	.....	19
Italiano	.....	25
Ελληνικά	.....	31
Deutsch	.....	37
Netherlands	.....	43
Interpretation	.....	49

# *Speed* Biogram™

## Veterinary Diagnostic kit

### For *in vitro* use only

#### ■ CLINICAL INTEREST

The causative organisms in urinary, dermatological or auricular infections are varied (bacteria, yeasts) and associations are possible. The blind use of antibiotic therapy is associated with the risk of therapeutic failure and may contribute to the selection of resistant strains.

The use of the Speed Biogram test on a simple liquid or cellular sample taken from the animal enables:

- the identification of pathogenic bacteria and/or yeast within 48h for dermatological, urinary, or auricular infections in dogs and cats.
- the determination of the antibiotic sensitivity profile of the organisms, using antibiotics that are available to the practitioner, within 24h to ensure the most appropriate treatment.

Speed Biogram takes into account the synergistic or antagonistic effects of various pathogenic agents, factors associated with the infected milieu (direct antibiotic sensitivity testing), as well as the concentration of organisms at the infected site (inoculum effect), to provide the most accurate simulation of *in vivo* conditions.

Moreover, the use of diagnostic tests and the visualisation of test results by the owner will help to ensure treatment compliance.

#### ■ PRINCIPLE

Speed Biogram is a diagnostic test for small animals; it is composed of a culture gallery comprising:

- 15 antibiotic wells enabling the determination of the sensitivity profile of the pathogens present in the sample.
- 6 wells for the identification of pathogenic bacteria.
- 1 well for the identification of *Malassezia* yeast.
- 2 control wells:
  - a ⊕ bacterial growth control well: if this well changes colour, presence of bacterial concentrations into the sample greater than 10<sup>3</sup> CFU/mL.
  - a ⊖ negative control well: if this well changes colour during the reading of the test, the test is invalid.

Speed Biogram is performed from a cellular (swab) or liquid sample. After inoculating the gallery, incubation at a controlled temperature (+37°C) enables the evaluation of the development of pathogenic microorganisms in the presence of antibiotics commonly used in veterinary practice and also the identification of the genus of the microorganism(s) present in the sample, via a simple colour change in each well.

#### ■ OPERATING PROCEDURE

##### ► FOR EACH TEST YOU WILL NEED:

One gallery, one swab, one bottle of *Preservative medium* (green cap), one bottle of *Culture medium* (Yellow cap), the *Staph supplement* bottle, the bottle of *Paraffin oil*,

one urine transfer pipette, one gallery holder, one incubator and a results sheet.  
**Never use reagents from different boxes.**

• **PRECAUTIONS:**

- It is advisable to wear protective gloves, protective clothing when taking the sample and during all stages of the test procedure (inoculation and reading).
- In case of contact of one reagent with skin, remove/take off immediately all contaminated clothing and rinse cautiously with water for several minutes.
- In case of contact of one reagent with eyes, rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- Perform the test procedures on a dry, clean, and disinfected work surface.

• **SAMPLING:**

Avoid applying antiseptic solution or local antibiotic before sampling. If you are taking a sample from an animal on local or systemic antibiotic treatment, it is advisable to take the sample *at least* 48 hours after the last administration.

*Cutaneous disease:*

- take the sample with a sterile swab.
- the sample should be taken from the walls of the lesions, after surgical debridement of any purulent accretions.

*Ear infections:*

- remove any superficial cerumen present on the walls of the auditory canal.
- take a deep swab of the walls of the canal using one of the sterile swabs provided.

*Urinary infections:*

- take a sterile urine sample by cystocentesis.

• **PREPARING THE SAMPLE:**

Sample on a swab:

- **Place the sample swab in the bottle of *Preservative medium*** (Green cap) and shake it vigorously in the bottle for a few seconds.
- Remove as much liquid as possible from the swab by pressing and turning the fibrous part against the walls of the bottle.
- Discard the swab in a biological waste container.
- Close the bottle of *Preservative medium* and shake well to homogenize the contents.

Urine sample:

- Use the urine transfer pipette to **place 2 drops of urine in the bottle of *Preservative medium*** (Green cap).
  - Discard the remainder of the sample in a biological waste container.
  - Close the bottle of *Preservative medium* and shake well to homogenize the contents.
- If you cannot perform the analysis immediately, the inoculated *Preservative medium* will remain stable for 48 hours at +4°C.

• **PREPARING THE GALLERY:**

- Using a gallery taken from its pouch just prior to use, record the name of the animal and the date of culture on the label. Remove the adhesive label covering the gallery.
- Stick the long edge of the label along the long edge of the gallery, to enable access

to all of the wells whilst keeping their identification visible.

- Please note that prior to bacteria growth the content of the control wells and of the antibiotic wells may not be visible with the naked eye.

• **INOCULATING THE CULTURE MEDIUM:**

- Through the included dropper cap, place **4 drops of the sowed Preservative medium** (Green cap) **into the bottle of Culture medium** (Yellow cap).

- Close the bottle of *Culture medium* and shake well to homogenize the contents.

• **INOCULATING THE GALLERY:**

- Through the included dropper cap, **distribute 3 drops of the inoculated Culture medium into each well** of the gallery.

- In the Staphylococcus identification well (STAPH well), add **2 drops of Staph supplement**.

- Then add **2 drops of Paraffin oil to each well, with the exception of the PSEUDO, E. COLI and STAPH wells**.

- Reposition the adhesive label on the gallery and make sure that you remove the transparent film from the back first for improved adherence.

- Place the gallery on the rigid card holder to improve the reading contrast and maintain an homogeneous temperature in all wells during culture.

• **CULTURE:**

After inoculation, immediately incubate the gallery at +37°C in an incubator.

• **READING:**

**1 - Reading the control wells:**

**Read control wells after 24 hours of incubation at +37°C. (plus/minus 2 hours)**

- The ⊖ negative control well should remain colourless:


- A change in the ⊕ growth control well from colourless to red, or with red streaks, indicates that the concentration of pathogenic bacteria is greater than 10<sup>3</sup> CFU/mL (no colour change of the well is observed in the presence of yeasts only)

Absence of pathogenic bacteria  
The ⊕ well remains colourless

Only yeast infection possible

Presence of pathogenic bacteria  
The ⊕ well turns red or contains red streaks

In this case, the antibiotic sensitivity profile can be read.

- The sample can be considered not infected by pathogenic bacteria if the ⊕ growth control well has not changed its color after 48 hours.

## 2 - Reading the antibiotic wells:

### Read antibiotic wells immediately after interpretation of the control wells

Antibiotic wells can be interpreted only if negative control well is colourless and red color is observed into positive control well.

No colour change.  
No bacterial growth.

Bacteria **SENSITIVE** to the antibiotic

or

Turns red or presence of red streaks.  
Bacterial growth

Bacteria **RESISTANT** to the antibiotic

## Antibiotic table of Speed Biogram:






























The identification of the antibiotic present in each well is printed on the sticky label, above each well.

Abbreviation	Antibiotic	Abbreviation	Antibiotic
AMO	Amoxicillin	SPI	Spiramycin
AMC	Amoxicillin + Clavulanic acid	CLI	Clindamycin
CFL	Cefalexin	NEO	Neomycin
CFT	Ceftiofur	GEN	Gentamicin
DOX	Doxycycline	SUL + TMP	Sulfonamides + Trimethoprim
FLU	Flumequine	FUS	Fusidic Acid
ENR	Enrofloxacin	PXB	Polymyxin B
MAR	Marbofloxacin		

## 3 – Reading the identification wells:

Read identification wells after 48 hours of incubation at +37°C, plus/minus 2 hours – (24h after control and antibiotic wells interpretation)

*Malassezia* yeast do not show any antibiotic profile on the gallery  
The test may reveal associations of different bacteria.

	Negative control	Positive control	ATB sensitive	ATB resistant	STAPH	STREP	PSEUDO	E. COLI	ENTERO BACT	PROTEUS	MALASSEZIA
<i>Staphylococcus spp.</i>	 Colourless	 Red	 Colourless	 Red	 Yellow	 Colourless	 Colourless	 Colourless	 Red	 Orange-Yellow	 Colourless
<i>Streptococcus spp.</i>	 Colourless	 Red	 Colourless	 Red	 Grey-blue	 Colourless	 Colourless	 Colourless	 Red	 Orange-Yellow	 Colourless
<i>Pseudomonas</i>	 Colourless	 Red	 Colourless	 Red	 Red	 Colourless	 Navy-blue	 Colourless	 Red	 Orange-Yellow	 Colourless
<i>E. coli</i>	 Colourless	 Red	 Colourless	 Red	 Red	 Colourless	 Colourless	 Blue-ring	 Orange-Yellow	 Orange-Yellow	 Colourless
<i>Proteus</i>	 Colourless	 Red	 Colourless	 Red	 Red	 Colourless	 Colourless	 Colourless	 Orange-Yellow	 Purple	 Colourless
<i>Enterobacteriaceae (except E.coli or Proteus)</i>	 Colourless	 Red	 Colourless	 Red	 Red	 Colourless	 Colourless	 Colourless	 Orange-Yellow	 Orange-Yellow	 Colourless
<i>Malassezia</i>	 Colourless	 Colourless	 Colourless	 Colourless	 Red	 Colourless	 Colourless	 Colourless	 Red	 Orange-Yellow	 Green-Brown

## Comments:

### Incomplete colour change:

If an identification well turns to an intermediate colour (often associated with turbid medium) after 48 hours of incubation at +37°C (e.g.: STAPH well turns from red to turbid orange), the well should be re-examined again 24h later (at 72 hours after inoculation)

### Result at 72 hours:

- Complete colour change: bacterial identification confirmed.
- Intermediate colour persists: bacterial identification is considered negative.

## ■ INTERPRETATION:

Results sheets are provided to record the identified bacteria and their antibiotic sensitivity and resistance profiles. The veterinarian's choice of antibiotic should take into account the results of the analysis and the pharmacokinetics and the toxicity of the chosen molecule with respect to the type and site of the infection.

## ■ RECOMMENDATIONS

### • STABILITY / STORAGE:

- The kit is stable at between +2°C and +8°C for 16 months from the date of manufacture (see expiry date on the label of the kit). Do not expose the kit to temperatures below freezing (0°C).
- It is advisable to leave all of the reagents and the gallery at room temperature for at least 15 minutes prior to use.
- Never use bottles of reagents from different boxes.

### • PRECAUTIONS OF USE:

- The volume and the number of drops of sample, *Preservative medium* or *Culture medium* recommended in the protocol are determining for the quality of the results.
- After sampling, microorganisms on a dry sample are stable for a maximum of 30 minutes. Immediately after sampling, empty the swab into a bottle of *Preservative medium*.
- Once inoculated, the *Preservative medium* can preserve the microorganisms present in the sample for:
  - ▶ 48 hours at +4°C
  - ▶ 4 to 6 months at -20°C

These recommendations are for guidance only; no diagnostic method is 100% accurate. The aim of this test is to identify the bacterial or fungal agent responsible for dermatological, auricular or urinary infections in dogs and cats, and to assist the veterinarian in their choice of treatment. The interpretation of the test by the veterinarian should always take into account the history, clinical examination, and the results of any further diagnostic tests. The definitive diagnosis is the prerogative and responsibility of the veterinarian.

Bio Veto Test and its distributors cannot be held responsible for any consequences linked to incorrect usage of this test or misinterpretation of the results.

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