




an- Test Strips

Ophthalmic Diagnostic Test Strips



Measurement of tear production
Each strip is marked with a measuring scale and impregnated with a blue or white indicator



Staining dead and degenerated cells and mucus.
Diagnose of KSC and superficial corneal epithelial abnormalities. Each strip is impregnated with 1.5 mg Rose Bengal



To detect corneal and conjunctival defects and to assess nasolacrimal patency.
Each strip is impregnated with Fluorescein Sodium



Staining of dead or degenerated conjunctival cells
Early detection of dry eye
Each strip is impregnated with 1.5 mg Lissamine Green